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Thank you for buying this Bosch appliance. These instructions contain valuable information about the unique features of this product.

Safety notes

Please read through the instructions carefully and proceed accordingly; store in a safe place!

This appliance is designed for domestic use or for use in household-like, non-commercial and/or medical environments.

Household-like environments include staff areas in shops, offices, agricultural and other small businesses or for use by guests at bed-and-breakfast establishments, small hotels and similar residential facilities.

Risk of slipping!

- Do not step on the scales with wet feet or if the surface of the scales is wet – Risk of slipping!
- Keep children under supervision to ensure that they do not play with the appliance.
- Do not let items fall on the scales.
- To prevent hazards, any repairs on this appliance must be made by our customer service.

These scales are to be used for the purpose of weighing people.

Overview

Illustrations **A** and **B**

- 1 Display
- 2 Battery compartment

Switching on the appliance

- Remove film from battery compartment **2** (Illustration **B**).
- Place the scales on a firm, dry and level surface.

Use

Illustration **C**

- To switch on the appliance, step on to the scales.
- Do not move until the display stops flashing. Distribute your weight evenly between both legs.
- The measurement can then be read. The scales show the weight in 100g increments.

Note: The scales will switch off automatically once the measurement has been taken.

Important: The scales automatically calibrate when switched off.

To ensure a correct reading, do not move the scales after switching off.

Cleaning the appliance

Never immerse the scales in water.

Do not use a steam cleaner.

Wipe the scales with a damp cloth; do not use scouring agents.

Error messages

If **Err** appears in the display, it is because:

- Too much weight (maximum: 180 kg).

If **0:00** appears, it is because:

- You moved during the measurement or did not stay on the scales long enough.
→ Stand still on the scales for at least five seconds.

If **Lo** appears, it is because:

- the batteries are too weak or dead.
→ Insert a fresh set of batteries (Illustration **B**).

Do not use rechargeable batteries! Use only batteries of the same type.

Disposal

Dispose of packaging in an environmentally-friendly manner. This appliance is labelled in accordance with European Directive 2012/19/EG concerning used electrical and electronic appliances (waste electrical and electronic equipment – WEEE). The guideline determines the framework for the return and recycling of used appliances as applicable throughout the EU. Please ask your specialist retailer about current disposal facilities.

Dispose of used batteries at public collection points or specialist dealers.

Guarantee

The guarantee conditions for this appliance are as defined by our representative in the country in which it is sold. Details regarding these conditions can be obtained from the dealer from whom the appliance was purchased. The bill of sale or receipt must be produced when making any claim under the terms of this guarantee.

Changes reserved.

Technical data

Dimensions (H x W x D)	25 x 300 x 300 mm
Batteries	CR2032
Measuring range in kg	kg (in steps of 0.1 kg) max. 180 kg

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GB Great Britain

BSH Home Appliances Ltd.
Grand Union House
Old Wolverton Road
Wolverton
Milton Keynes MK12 5PT
To arrange an engineer visit, to
order spare parts and accessories
or for product advice please visit
www.bosch-home.co.uk
or call
Tel.: 0844 892 8979*

*Calls from a BT landline will be charged at
up to 3 pence per minute. A call set-up fee
of up to 6 pence may apply.

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IE Republic of Ireland

BSH Home Appliances Ltd.
Unit F4, Ballymount Drive
Ballymount Industrial Estate
Walkinstown
Dublin 12
Service Requests,
Spares and Accessories
Tel.: 01450 2655
Fax: 01450 2520
www.bosch-home.co.uk

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IN India, Bhārat, भारत

BSH Customer Service Front Office
Shop No.4, Everest Grande,
Opp. Shanti Nagar Bus Stop,
Mahakali Caves Road, Andheri East
Mumbai 400093

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Robert Bosch Hausgeräte GmbH