## PIL1 /PIL2

de Gebrauchsanleitung  
en Operating instructions  
fr Notice d'utilisation  
it Istruzioni per l'uso  
nl Gebruiksaanwijzing  
da Brugsanvisning  
no Bruksanvisning  
sv Bruksanvisning  
fi Käyttöohje  
es Instrucciones de uso

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**Pay attention!**  
**Read operating instructions**  
**and keep in a safe place!**

**Application part of type B**

**Protection class II**

**Alternating current**

## General safety instructions

This infrared radiator is intended for domestic use only and must not be used for industrial or medical applications.

**Connect and operate the appliance only in accordance with the specifications on the rating plate.**

**This infrared radiator satisfies the CE regulations and complies with the directives 89/336/EG and 93/42/EG.**

### **⚠ Risk of electric shock!**

Do not use the appliance if the power cord and/or appliance are damaged.

After the use and if a fault occurs, immediately pull out the mains plug.

To protect the user, all repairs to the appliance, e.g. replacement of the power cord, should be carried out only by our customer service.

Do NOT

- bring the power cord into contact with hot parts,
- pull the power cord over sharp edges,
- carry the appliance by the power cord.

Never bring the infrared radiator into contact with water, other liquids or humidity.

If the infrared radiator falls under or into water, stop using the appliance and have it checked by an authorised service centre.

To ensure that the power supply has been disconnected, pull the mains plug out of the socket.

Do not place the infrared radiator on hot surfaces or near heat sources.

### **⚠ Risk of overheating!**

Do not leave the infrared radiator unattended while on.

Do not place the infrared radiator near highly flammable objects (**minimum distance 80 cm**) and never cover.

Do not hang the infrared radiator from the wall or ceiling.

### **⚠ Risk of injury!**

Children or people whose physical, sensory or mental abilities prevent them from using the infrared radiator safely must **not use the infrared radiator without help or supervision.**

When treating your face, always close your eyes; never look directly into the light emitted by the infrared bulb.

Do not touch the infrared radiator while in use (becomes hot). First, leave it to cool down completely.

## Special safety instructions

**⚠ In the case of actual illness or skin disease or if harmless looking symptoms persist, please consult a doctor.**

**⚠ Do not use the infrared radiator if taking pain-relieving medication or after drinking alcohol as this could restrict thermal sensitivity.**  
**Avoid cooling the treated part of the body too quickly.**

Use the infrared radiator with an infrared bulb only.

Do not subject the infrared bulb to violent impacts (e.g. by banging down the infrared radiator).

## Application options

The skin is heated by irradiation with infrared light. This significantly increases the blood supply to the skin and stimulates metabolism. In particular, the infrared radiator aids healing processes and alleviates pain.

Infrared light is an established means of preventing and treating colds and alleviating muscular pain. It is also a facial and beauty care aid and a particularly effective treatment for unhealthy skin.

## Use

- Select a distance between the infrared bulb and your body to ensure that an agreeable and pleasant heat is irradiated.

Recommendation:

- Distance approx. 40–60 cm
- Treatment duration max. 15 minutes
- 4 Thanks to 4 tilt positions the infrared radiator can be set at the required position to suit the individual (Fig. 1).
- After use, always pull out the mains plug.

### Additional information concerning model PIL 2200

#### Fig. 2

This infrared radiator also features a timer.

Control panel:

- 0 Off
- 1 Continuous operation
- 2 Irradiation duration max. 15 minutes

## Changing the infrared bulb

- Pull out the mains plug, leave the infrared radiator to cool down.
- Unscrew the infrared bulb and fit an appropriate replacement bulb.

Replacement bulb: use original type only!

- PIL 1000: R95E/100 W
- PIL 2200: PAR38E/150 W

## Cleaning, storage

### Caution:

**Before cleaning, always pull out the mains plug and leave the infrared radiator to cool down!**

- Clean the infrared radiator with a damp cloth.  
Do not use detergents or solvents. These may damage the appliance.  
Do not use a steam cleaner!
- Wait until the infrared radiator has cooled down before putting it away.
- On model PIL 2200 the power cord can be inserted into a cord store (Fig. 3).

## Disposal

Please ask your dealer or inquire at your local authority about current means of disposal.

## Guarantee

The guarantee conditions for this appliance are as defined by our representative in the country in which it is sold. Details regarding these conditions can be obtained from the dealer from whom the appliance was purchased.

The bill of sale or receipt must be produced when making any claim under the terms of this guarantee.

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**Robert Bosch Hausgeräte GmbH**
Luxoplast GmbH