

### Long-life recirculation kit, initial equipment DSZ1WW1J1



## Long Life Regenerative filter: Removes odors reliably over a period of up to 10 years.

- For recirculation operation
- Fits 80cm Glassdraft
- Can be used multiple times, the regenerative active carbon filter is regenerated in the oven with a recommendation of every 4 months. Simply place in the oven for 2 hours at 200°C.
- The Long Life Regenerative recirculation kit includes 2 filter cassettes, a flexible aluminium ducting and a ducting clamp

### **Technical Data**

270x350x450 mm
180.0 x 80.0 x 120.0
24
4.1 kg
5.2 kg
4242005241644





### Long-life recirculation kit, initial equipment DSZ1WW1J1

# Long Life Regenerative filter: Removes odors reliably over a period of up to 10 years.

### **Optional accessories**

- For recirculation operation
- Fits 80cm Glassdraft
- Can be used multiple times, the regenerative active carbon filter is regenerated in the oven with a recommendation of every 4 months. Simply place in the oven for 2 hours at 200°C.
- The Long Life Regenerative recirculation kit includes 2 filter cassettes, a flexible aluminium ducting and a ducting clamp